FundusAnalysis.py

```
class FundusAnalysis:
    def initial(self, record_dir, model_dir):
        """
        Initial models and config. A time-consuming process.
        Args:
            record_dir: str. The path to store output results.
            model_dir: str. The path to load models.
        """

    def analysis(self, images, va=True, ignore_exception):
        """
        Analyse batch of images.
        Args:
            images: list(str). The path of images to analysis.
            va: bool. If do the vessel analysis.
            ignore_exception: bool. If continue the batch process when any exception
            raises. True for application and False for testing.
        Returns:
            reports: list(dict). List of report. The detail of key and value is shown
            below.
            exceptions: dict. Key is image_path and value is exception traceback if any
            exception raises.
            statuses: dict. Key is image_path and the value is status code.
            Status Code:
                0: Success
                1: Not Fundus Images
                2: Optical Disc Unclear
        """

    def _analysis_one_shot(self, va, image_path, image_name):
        """
        Private function. Should not be called.
        Args:
            va: bool. If do the vessel analysis.
            image_path: str. The path of image.
            image_name: str. The name of image shown in the log.
        Returns:
            report: dict. The result of analysis.
        """
```

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